

Product Datasheet

Human Protachykinin-1 ELISA Kit (orb1213796)

Description	Human Protachykinin-1, TAC1 ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	TAC1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	15,003
Uniprot ID	P20366
Solubility (25°C)	6400ng/L
Sensitivity	10.47ng/L
Entrez	6863
Expiration Date	Please enquire.