

Product Datasheet

Human Fibrillin-1 ELISA Kit (orb1213805)

Description	Human Fibrillin-1, FBN1 ELISA Kit
Reactivity	Human
Range	0.5-150ng/ml
Tested Applications	ELISA
Target	FBLN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	77,214
Uniprot ID	P23142
Solubility (25°C)	160ng/ml
Sensitivity	0.18ng/ml
Entrez	2192
Expiration Date	Please enquire.