

Product Datasheet

Human Selenoprotein-S ELISA Kit (orb1213806)

Description	Human Selenoprotein-S, SE-S ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	SELENOS
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,163
Uniprot ID	Q9BQE4
Solubility (25°C)	4800ng/L
Sensitivity	10.39ng/L
Entrez	55829
Expiration Date	Please enquire.