

Product Datasheet

Human Annexin A1 ELISA Kit (orb1213810)

Description	Human Annexin A1, ANXA1 ELISA Kit
Reactivity	Human
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	ANXA1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,714
Uniprot ID	P04083
Solubility (25°C)	96ng/ml
Sensitivity	0.18ng/ml
Entrez	301
Expiration Date	Please enquire.