

Product Datasheet

Human Retinal Dehydrogenase 1 ELISA Kit (orb1213823)

Description	Human Retinal Dehydrogenase 1, ALDH1A1 ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	ALDH1A1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,862
Uniprot ID	P00352
Solubility (25°C)	80ng/ml
Sensitivity	0.11ng/ml
Entrez	216
Expiration Date	Please enquire.