

Product Datasheet

Human Selenium-binding Protein 1 ELISA Kit (orb1213830)

Description	Human Selenium-binding Protein 1, SELENBP1 ELISA Kit
Reactivity	Human
Range	70-15000ng/L
Tested Applications	ELISA
Target	CABP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,838
Uniprot ID	Q9NZU7
Solubility (25°C)	16000ng/L
Sensitivity	32.95ng/L
Entrez	9478
Expiration Date	Please enquire.