

Product Datasheet

Human Nuclear Factor Erythroid 2-related Factor 2 ELISA Kit (orb1213832)

Description	Human Nuclear Factor Erythroid 2-related Factor 2, NRF2 ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	NFE2L2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	67,827
Uniprot ID	Q16236
Solubility (25°C)	80ng/ml
Sensitivity	0.11ng/ml
Entrez	4780
Expiration Date	Please enquire.