

Product Datasheet

Human Cyclic Amp-dependent Transcription Factor Atf-2 ELISA Kit (orb1213834)

Description	Human Cyclic Amp-dependent Transcription Factor Atf-2, ATF2 ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	ATF2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,537
Uniprot ID	P15336
Solubility (25°C)	16ng/ml
Sensitivity	0.024ng/ml
Entrez	1386
Expiration Date	Please enquire.