

Product Datasheet

Human Fructosamine-3-kinase ELISA Kit (orb1213837)

Description	Human Fructosamine-3-kinase, FN3K ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	FN3K
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,171
Uniprot ID	Q9H479
Solubility (25°C)	48ng/ml
Sensitivity	0.054ng/ml
Entrez	64122
Expiration Date	Please enquire.