

Product Datasheet

Human Cytochrome P450 1A1 ELISA Kit (orb1213841)

Description	Human Cytochrome P450 1A1, CYP1A1 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	CYP1A1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	58,165
Uniprot ID	P04798
Solubility (25°C)	2400ng/L
Sensitivity	5.25ng/L
Entrez	1543
Expiration Date	Please enquire.