

Product Datasheet

Human Ninein/Gsk3B Interacting Protein ELISA Kit (orb1213849)

Description	Human Ninein/Gsk3B Interacting Protein, NIN ELISA Kit
Reactivity	Human
Range	0.1-30ng/ml
Tested Applications	ELISA
Target	GSKIP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	15,648
Uniprot ID	Q9P0R6
Solubility (25°C)	32ng/ml
Sensitivity	0.05ng/ml
Entrez	51527
Expiration Date	Please enquire.