

Product Datasheet

Human Formin-binding Protein 1-like ELISA Kit (orb1213854)

Description	Human Formin-binding Protein 1-like, FBNP1L ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	FBNP1L
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	70,065
Uniprot ID	Q5TON5
Solubility (25°C)	3200ng/L
Sensitivity	7.98ng/L
Entrez	54874
Expiration Date	Please enquire.