

Product Datasheet

Human Spry Domain-containing Socs Box Protein 1 ELISA Kit (orb1213860)

Description	Human Spry Domain-containing Socs Box Protein 1, SPSB1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	SPSB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,942
Uniprot ID	Q96BD6
Solubility (25°C)	4800ng/L
Sensitivity	9.96ng/L
Entrez	80176
Expiration Date	Please enquire.