

Product Datasheet

Human Cytokine Receptor-like Factor 2 ELISA Kit (orb1213861)

Description	Human Cytokine Receptor-like Factor 2, CRLF2 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	CRLF2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,013
Uniprot ID	Q9HC73
Solubility (25°C)	24ng/ml
Sensitivity	0.026ng/ml
Entrez	64109
Expiration Date	Please enquire.