

Product Datasheet

Human Atrial Natriuretic Peptide Receptor 2 ELISA Kit (orb1213863)

Description	Human Atrial Natriuretic Peptide Receptor 2, NPR2 ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	NPR2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	117,022
Uniprot ID	P20594
Solubility (25°C)	12.8ng/ml
Sensitivity	0.014ng/ml
Entrez	4882
Expiration Date	Please enquire.