

Product Datasheet

Human Cytochrome P450 3A4 ELISA Kit (orb1213896)

Description	Human Cytochrome P450 3A4, CYP3A4 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	CYP3A4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	57,343
Uniprot ID	P08684
Solubility (25°C)	3200ng/L
Sensitivity	6.68ng/L
Entrez	1576
Expiration Date	Please enquire.