

Product Datasheet

Human Zinc-alpha-2-glycoprotein ELISA Kit (orb1213905)

Description	Human Zinc-alpha-2-glycoprotein, ZAG ELISA Kit
Reactivity	Human
Range	3-900ug/ml
Tested Applications	ELISA
Target	AZGP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	34,259
Uniprot ID	P25311
Solubility (25°C)	960µg/ml
Sensitivity	1.49µg/ml
Entrez	563
Expiration Date	Please enquire.