

Product Datasheet

Human Annexin A5 ELISA Kit (orb1213926)

Description	Human Annexin A5, ANXA5 ELISA Kit
Reactivity	Human
Range	0.1-35ng/ml
Tested Applications	ELISA
Target	ANXA5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,937
Uniprot ID	P08758
Solubility (25°C)	40ng/ml
Sensitivity	0.049ng/ml
Entrez	308
Expiration Date	Please enquire.