

Product Datasheet

Human Suppressors of Cytokine Signaling 1 ELISA Kit (orb1213933)

Description	Human Suppressors of Cytokine Signaling 1, SOCS1 ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	SOCS1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,551
Uniprot ID	015524
Solubility (25°C)	80ng/ml
Sensitivity	0.102ng/ml
Entrez	8651
Expiration Date	Please enquire.