

Product Datasheet

Human Calreticulin 3 ELISA Kit (orb1213939)

Description	Human Calreticulin 3, CALR3 ELISA Kit
Reactivity	Human
Range	0.5-150ng/ml
Tested Applications	ELISA
Target	CALR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	48,142
Uniprot ID	P27797
Solubility (25°C)	160ng/ml
Sensitivity	0.24ng/ml
Entrez	811
Expiration Date	Please enquire.