

Product Datasheet

Human Peroxiredoxin 2 ELISA Kit (orb1213950)

Description	Human Peroxiredoxin 2, PRDX2 ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	PRDX2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,892
Uniprot ID	P32119
Solubility (25°C)	6400ng/L
Sensitivity	9.58ng/L
Entrez	7001
Expiration Date	Please enquire.