

Product Datasheet

Human Calbindin ELISA Kit (orb1213976)

Description	Human Calbindin, CALB ELISA Kit
Reactivity	Human
Range	3-900ng/L
Tested Applications	ELISA
Target	CALB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,025
Uniprot ID	P05937
Solubility (25°C)	960ng/L
Sensitivity	1.49ng/L
Entrez	793
Expiration Date	Please enquire.