

Product Datasheet

Human Myc Box-dependent-interacting Protein 1 ELISA Kit (orb1213981)

Description	Human Myc Box-dependent-interacting Protein 1, BIN1 ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	BIN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	64,699
Uniprot ID	000499
Solubility (25°C)	4000ng/L
Sensitivity	9.97ng/L
Entrez	274
Expiration Date	Please enquire.