

Product Datasheet

Human Breast Cancer Anti Estrogen Resistance 1 ELISA Kit (orb1213982)

Description	Human Breast Cancer Anti Estrogen Resistance 1, BCAR1 ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	BCAR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	93,372
Uniprot ID	P56945
Solubility (25°C)	6400ng/L
Sensitivity	10.26ng/L
Entrez	9564
Expiration Date	Please enquire.