

Product Datasheet

Human Biliverdin Reductase A ELISA Kit (orb1213983)

Description	Human Biliverdin Reductase A, BLVRA ELISA Kit
Reactivity	Human
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	BLVRA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	33,428
Uniprot ID	P53004
Solubility (25°C)	16ng/ml
Sensitivity	0.023ng/ml
Entrez	644
Expiration Date	Please enquire.