

Product Datasheet

Human Battenin ELISA Kit (orb1213984)

Description	Human Battenin, BSND ELISA Kit
Reactivity	Human
Range	20-3800ng/L
Tested Applications	ELISA
Target	CLN3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	47,623
Uniprot ID	Q13286
Solubility (25°C)	4000ng/L
Sensitivity	9.64ng/L
Entrez	1201
Expiration Date	Please enquire.