

Product Datasheet

Human Protein Atonal Homolog 1 ELISA Kit (orb1213985)

Description	Human Protein Atonal Homolog 1, ATOH1 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	ATOH1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,160
Uniprot ID	Q92858
Solubility (25°C)	2400ng/L
Sensitivity	4.96ng/L
Entrez	474
Expiration Date	Please enquire.