

## Product Datasheet

### Human Aryl Hydrocarbon Receptor Repressor ELISA Kit (orb1213989)

|                                |                                                                 |
|--------------------------------|-----------------------------------------------------------------|
| <b>Description</b>             | Human Aryl Hydrocarbon<br>Receptor Repressor, AHRR<br>ELISA Kit |
| <b>Reactivity</b>              | Human                                                           |
| <b>Range</b>                   | 10-2000ng/L                                                     |
| <b>Tested<br/>Applications</b> | ELISA                                                           |
| <b>Target</b>                  | AHRR                                                            |
| <b>Storage</b>                 | 2-8°C                                                           |
| <b>Note</b>                    | For research use only                                           |
| <b>Application notes</b>       | Sandwich                                                        |
| <b>Sample Types</b>            | Serum, plasma, cell culture<br>supernates                       |
| <b>MW</b>                      | 76,265                                                          |
| <b>Uniprot ID</b>              | <b>A9YTQ3</b>                                                   |
| <b>Solubility (25°C)</b>       | 2400ng/L                                                        |
| <b>Sensitivity</b>             | 5.35ng/L                                                        |
| <b>Entrez</b>                  | <b>57491</b>                                                    |
| <b>Expiration Date</b>         | Please enquire.                                                 |