

Product Datasheet

Human Adipocyte Enhancer-binding Protein 1 ELISA Kit (orb1213994)

Description	Human Adipocyte Enhancer-binding Protein 1, AEBP1 ELISA Kit
Reactivity	Human
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	AEBP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	130,929
Uniprot ID	Q8IUX7
Solubility (25°C)	12.8ng/ml
Sensitivity	0.02ng/ml
Entrez	165
Expiration Date	Please enquire.