

Product Datasheet

Human Acyloxyacyl Hydrolase ELISA Kit (orb1213995)

Description	Human Acyloxyacyl Hydrolase, AOA ELISA Kit
Reactivity	Human
Range	0.05-20ng/L
Tested Applications	ELISA
Target	AOAH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	65,105
Uniprot ID	P28039
Solubility (25°C)	24ng/L
Sensitivity	0.025ng/L
Entrez	313
Expiration Date	Please enquire.