

Product Datasheet

Human Cystatin-Sn ELISA Kit (orb1213998)

Description	Human Cystatin-Sn, CST1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	CST1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,388
Uniprot ID	P01037
Solubility (25°C)	24ng/ml
Sensitivity	0.023ng/ml
Entrez	1469
Expiration Date	Please enquire.