

Product Datasheet

Human Inter-alpha-trypsin Inhibitor Heavy Chain H1 ELISA Kit (orb1214013)

Description	Human Inter-alpha-trypsin Inhibitor Heavy Chain H1, ITIH1 ELISA Kit
Reactivity	Human
Range	10-2000ng/ml
Tested Applications	ELISA
Target	ITIH1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	101,389
Uniprot ID	P19827
Solubility (25°C)	2400ng/ml
Sensitivity	4.98ng/ml
Entrez	3697
Expiration Date	Please enquire.