

Product Datasheet

Human Golgi Phosphoprotein 3 ELISA Kit (orb1214017)

Description	Human Golgi Phosphoprotein 3, GOLPH3 ELISA Kit
Reactivity	Human
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	GOLPH3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	33,811
Uniprot ID	Q9H4A6
Solubility (25°C)	64ng/ml
Sensitivity	0.11ng/ml
Entrez	64083
Expiration Date	Please enquire.