

Product Datasheet

Human Glycogenin 1 ELISA Kit (orb1214020)

Description	Human Glycogenin 1, GYG1 ELISA Kit
Reactivity	Human
Range	0.5-100ng/ml
Tested Applications	ELISA
Target	GYG1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,384
Uniprot ID	P46976
Solubility (25°C)	128ng/ml
Sensitivity	0.28ng/ml
Entrez	2992
Expiration Date	Please enquire.