

Product Datasheet

Human Prosaposin Receptor Gpr37 ELISA Kit (orb1214027)

Description	Human Prosaposin Receptor Gpr37, GPR37 ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	GPR37
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	67,114
Uniprot ID	O15354
Solubility (25°C)	32ng/ml
Sensitivity	0.027ng/ml
Entrez	2861
Expiration Date	Please enquire.