

Product Datasheet

Human Fibroblast Growth Factor 3 ELISA Kit (orb1214032)

Description	Human Fibroblast Growth Factor 3, FGF3 ELISA Kit
Reactivity	Human
Range	5-1500ng/L
Tested Applications	ELISA
Target	FGF3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,887
Uniprot ID	P11487
Solubility (25°C)	1600ng/L
Sensitivity	2.64ng/L
Entrez	2248
Expiration Date	Please enquire.