

## Product Datasheet

### Human Disks Large-associated Protein 5 ELISA Kit (orb1214040)

|                            |                                                          |
|----------------------------|----------------------------------------------------------|
| <b>Description</b>         | Human Disks Large-associated Protein 5, DLGAP5 ELISA Kit |
| <b>Reactivity</b>          | Human                                                    |
| <b>Range</b>               | 0.5-300ng/ml                                             |
| <b>Tested Applications</b> | ELISA                                                    |
| <b>Target</b>              | DLGAP5                                                   |
| <b>Storage</b>             | 2-8°C                                                    |
| <b>Note</b>                | For research use only                                    |
| <b>Application notes</b>   | Sandwich                                                 |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                   |
| <b>MW</b>                  | 95,115                                                   |
| <b>Uniprot ID</b>          | <b>Q15398</b>                                            |
| <b>Solubility (25°C)</b>   | 320ng/ml                                                 |
| <b>Sensitivity</b>         | 0.31ng/ml                                                |
| <b>Entrez</b>              | <b>9787</b>                                              |
| <b>Expiration Date</b>     | Please enquire.                                          |