

Product Datasheet

Human Protein Diaphanous Homolog 1 ELISA Kit (orb1214043)

Description	Human Protein Diaphanous Homolog 1, DIAPH1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	DIAPH1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	141,347
Uniprot ID	O60610
Solubility (25°C)	24ng/ml
Sensitivity	0.019ng/ml
Entrez	1729
Expiration Date	Please enquire.