

Product Datasheet

Human Desmocollin 1 ELISA Kit (orb1214048)

Description	Human Desmocollin 1, DSC1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	DSC1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	99,987
Uniprot ID	Q08554
Solubility (25°C)	3200ng/L
Sensitivity	6.55ng/L
Entrez	1823
Expiration Date	Please enquire.