

Product Datasheet

Human Torsin 1A ELISA Kit (orb1214061)

Description	Human Torsin 1A, TOR1A ELISA Kit
Reactivity	Human
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	TOR1A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	37,809
Uniprot ID	O14656
Solubility (25°C)	32ng/ml
Sensitivity	0.022ng/ml
Entrez	1861
Expiration Date	Please enquire.