

Product Datasheet

Human Tnf Receptor Associated Factor 6 ELISA Kit (orb1214062)

Description	Human Tnf Receptor Associated Factor 6, TRAF6 ELISA Kit
Reactivity	Human
Range	20-3500ng/L
Tested Applications	ELISA
Target	TRAF6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	59,573
Uniprot ID	Q9Y4K3
Solubility (25°C)	4000ng/L
Sensitivity	9.57ng/L
Entrez	7189
Expiration Date	Please enquire.