

Product Datasheet

Human Statherin ELISA Kit (orb1214070)

Description	Human Statherin, STATH ELISA Kit
Reactivity	Human
Range	30-9000ng/ml
Tested Applications	ELISA
Target	STATH
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	7,304
Uniprot ID	P02808
Solubility (25°C)	9600ng/ml
Sensitivity	16.42ng/ml
Entrez	6779
Expiration Date	Please enquire.