

Product Datasheet

Human Semaphorin 6D ELISA Kit (orb1214077)

Description	Human Semaphorin 6D, SEMA6D ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	SEMA6D
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	119,872
Uniprot ID	Q8NFY4
Solubility (25°C)	2400ng/L
Sensitivity	4.45ng/L
Entrez	80031
Expiration Date	Please enquire.