

Product Datasheet

Human Coronin 1C ELISA Kit (orb1214093)

Description	Human Coronin 1C, CORO1C ELISA Kit
Reactivity	Human
Range	20-6000ng/L
Tested Applications	ELISA
Target	CORO1C
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,249
Uniprot ID	Q9ULV4
Solubility (25°C)	6400ng/L
Sensitivity	9.19ng/L
Entrez	23603
Expiration Date	Please enquire.