

Product Datasheet

Human Coronin 1A ELISA Kit (orb1214095)

Description	Human Coronin 1A, CORO1A ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	CORO1A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,026
Uniprot ID	P31146
Solubility (25°C)	24ng/ml
Sensitivity	0.026ng/ml
Entrez	11151
Expiration Date	Please enquire.