

Product Datasheet

Human Proteasome Subunit Beta Type 6 ELISA Kit (orb1214103)

Description	Human Proteasome Subunit Beta Type 6, PSMB6 ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	PSMB6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,358
Uniprot ID	P28072
Solubility (25°C)	1600ng/L
Sensitivity	2.98ng/L
Entrez	5694
Expiration Date	Please enquire.