

Product Datasheet

Human Neogenin 1 ELISA Kit (orb1214128)

Description	Human Neogenin 1, NEO1 ELISA Kit
Reactivity	Human
Range	5-1000ng/L
Tested Applications	ELISA
Target	NEUROG1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,718
Uniprot ID	Q92886
Solubility (25°C)	1280ng/L
Sensitivity	3.11ng/L
Entrez	4762
Expiration Date	Please enquire.