

Product Datasheet

Human Protein Wnt-2 ELISA Kit (orb1214136)

Description	Human Protein Wnt-2, WNT2 ELISA Kit
Reactivity	Human
Range	0.2-60ng/mL
Tested Applications	ELISA
Target	WNT2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,418
Uniprot ID	P09544
Solubility (25°C)	64ng/mL
Sensitivity	0.091ng/mL
Entrez	7472
Expiration Date	Please enquire.