

Product Datasheet

Human Protein Wnt-3A ELISA Kit (orb1214137)

Description	Human Protein Wnt-3A, WNT3A ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	WNT3A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,365
Uniprot ID	P56704
Solubility (25°C)	24ng/ml
Sensitivity	0.02ng/ml
Entrez	89780
Expiration Date	Please enquire.