

Product Datasheet

Human Serine ELISA Kit (orb1214139)

Description	Human Serine, Threonine-protein Kinase Pim-1, PIM1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	PIM1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,686
Uniprot ID	P11309
Solubility (25°C)	24ng/ml
Sensitivity	0.025ng/ml
Entrez	5292
Expiration Date	Please enquire.